

Siemens Healthcare Diagnostics GmbH, SHS EMEA CEET QT, Siemensstrasse 90,
1210 Vienna

Name	M.A. Roland Ertl
Department	SHS EMEA CEET QT
Telephone	+43 51707-38274
Mobile	+43 (664) 8011738274
E-mail	roland.re.ertl@siemens-healthineers.com
Date	December 23, 2020
Document Ref#	PH-21-001 A OUS

Urgent Field Safety Notice:

Atellica® COAG 360 System Software Version 1.8

Automatic re-dilutions not measured

Dear Sirs,

Our records indicate that your facility may have received the following product:

Table 1. Atellica COAG 360 Systems Affected Product(s)

Atellica COAG 360 System Software Version	Siemens Material Number (SMN)
1.8	11468626

Reason for Correction

The purpose of this communication is to inform you of an issue with the product indicated in Table 1 above and provide instructions on actions that your laboratory must take.

Siemens Healthcare Diagnostics Products GmbH has confirmed that the Atellica COAG 360 Software Version 1.8 may cause a delay in patient sample result reporting as the automatic re-dilutions are not measured under certain conditions:

- A sample result fulfills the criterion for a required automatic re-dilution
- Request for the automatic re-dilution is set and no further measurement is placed in the Joblist or in the Control Journal

Siemens Healthcare Diagnostics Products GmbH investigated the issue and determined as root cause a software bug and the solution will be provided with the upcoming Atellica COAG 360 System Software Version 1.9.

Please see section Additional Information, which provides a workaround to avoid not performed automatic re-dilutions.

Risk to Health

A low health risk is associated to this issue as due to the lack of information the above described issue may lead to a delay in reporting of patient sample results for emergency / urgency samples.

As no results were obtained in the failure mode, it is not necessary to review previous test results or to repeat testing. Therefore, also no look back is required.

Actions to be Taken by the Customer

- Please review this letter with your Medical Director.
- Perform the instructions provided in Additional Information.
- Complete and return the Field Correction Effectiveness Check Form attached to this letter within 30 days.
- Please retain this letter with your laboratory records and forward this letter to those who may have received this product.

Additional Information

Following procedure describes the workaround to avoid not performed automatic re-dilutions (see Table 2, which lists the methods which may initiate an automatic re-dilution if required and are therefore influenced by the issue).

- If the automatic re-dilution is the last measurement order, this re-dilution will not be started.
- Select any other measurement request in the Joblist / Control Journal or place a new sample / control vial (with open orders) in the Sample Manager.
- As soon as the Atellica COAG 360 System starts to plan this ordered measurement, the open re-dilution testing will be performed as well, and a corresponding finalization time will be displayed.
- To save sample material and reagent, the additional request that was used to trigger to start the measurement of the automatic re-dilution measurement may be deleted afterwards, if this was not required.
- Please Note: A deletion of the remaining automatic re-dilution request is not recommended, as the initial testing will be repeated, and the issue of the delayed automatic re-dilution might happen again.

Table 2. Atellica COAG 360 Systems Affected Method(s)

Assay Name	Assay No.	Condition for re-Dilution	Alternative Dilution	
Antithrombin Ag	14200	If result is below calibration range use dilution	1:	2
		If result is above calibration range use dilution	1:	8
Antithrombin Berichrom (Anti IIa)	14000	If result is above calibration range use dilution	1:	12
C1-Inhibitor Berichrom	33100	If result is below calibration range use dilution	1:	2
		If result is above calibration range use dilution	1:	8
D-Dimer INNOVANCE mg/L FEU	15000	If result is above calibration range use dilution	1:	10
D-Dimer INNOVANCE µg/L FEU	15010	If result is above calibration range use dilution	1:	10
D-Dimer LOCI µg/L FEU	71000	If result is above 9000 µg/L FEU use dilution	1:	3
D-Dimer LOCI mg/L FEU	71001	If result is above 9 mg/L FEU use dilution	1:	3
Factor IX Actin FS	27100	If result is above calibration range use dilution	1:	16
Factor IX Actin FS 1:2	27101	If result is above calibration range use dilution	1:	16
Factor IX Actin FS 1:4	27102	If result is above calibration range use dilution	1:	16
Factor IX Actin FSL	27200	If result is above calibration range use dilution	1:	40
Factor IX Pathromtin SL	27000	If result is above calibration range use dilution	1:	40
Factor VIII Actin FS	26100	If result is above calibration range use dilution	1:	16
Factor VIII Actin FS 1:2	26100	If result is above calibration range use dilution	1:	16
Factor VIII Actin FS 1:4	26100	If result is above calibration range use dilution	1:	16
Factor VIII Actin FSL	26200	If result is above calibration range use dilution	1:	40
Factor VIII Pathromtin SL	2600	If result is above calibration range use dilution	1:	40
Factor XI Actin FS	29100	If result is above calibration range use dilution	1:	20
Fibrinogen Multifibren U	13100	If result is above calibration range use dilution	1:	2
Fibrinogen Thrombin Reagent	13000	If result is below 1,05 g/L use dilution	1:	1
		If result is above 4,8 g/L use dilution	1:	6
FVIII chromogenic	26400	If result is above calibration range use dilution	1:	30
Von WillebrandFactor (vWF Ac)	16001	If result is above 150 % use dilution	1:	12
Von WillebrandFactor (vWF Ag)	16000	If result is above 170 % use dilution	1:	8

We apologize for the inconvenience this situation may cause. If you have any questions, please contact your Siemens Healthineers Remote Services Center of your local Siemens Healthineers Technical Support Representative.

Letter of December 23, 2020
to



Sincerely yours,

Siemens Healthcare Diagnostics GmbH

A handwritten signature in black ink, appearing to read "Franz Schwarz".

i.V. Dipl. Ing. Franz Schwarz
Quality Management CEE

A handwritten signature in black ink, appearing to read "Gottfried Pechtl".

i.A. DI Gottfried Pechtl
Product Manager Austria & SEE

—