# IMPORTANT RECALL NOTICE



Innovation. Experience. Performance.



DATE: 04/23/2018

**Dear Valued Customers:** 

<u>Chart has received reports of a vacuum leak or failure that could compromise the product and may affect Aluminum Dewar models listed below.</u>

The reliability, safety, and quality of our products is of utmost importance to Chart Industries. In addition to this communication, we will personally call you to verify receipt of this recall notice. We take seriously the trust you and your customers place in the quality of Chart's products. Please contact us at <a href="MVEtankservice@chartindustries.com">MVEtankservice@chartindustries.com</a> if you have any questions or comments. Thank you for your business and cooperation.

#### **SERIAL NUMBERS AFFECTED:**

Aluminum Dewars with Serial Number Range: NPB201750xxxxx through NPB201752xxxxx and NPB201801xxxxx through NPB201810xxxxx (manufactured Week 50 of 2017 through Week 10 of 2018).

A list further describing potentially affected Model numbers and names is provided at the end of this recall notice.

This recall does not pertain to Chart products manufactured anywhere other than New Prague, Minnesota, USA. This recall does not pertain to Chart's Stainless Steel Cryogenic freezers.

#### **ISSUE:**

Chart is presently investigating the possible cause of the VACUUM LEAK AND/OR FAILURE which may be due to inadequate adhesion of the composite neck to the aluminum unit. Based on information to date, this issue appears to be an isolated occurrence involving the machine and binding agent used during the manufacturing process. Chart is recalling with immediate effect the affected models and will replace affected units.

### **REQUESTED IMMEDIATE ACTION:**

If you have a unit in question (within the serial number range noted above) you must:

• <u>Immediately remove the unit(s) from service;</u>

AND, EITHER:

- Return to Chart Inc.; OR
- After following the standard procedures to empty the unit, drill a hole in the retired unit near the serial number. Drilling the unit will assure the unit is not returned to service. Email or transmit photographic evidence to Chart of compliance with this Recall notice, depicting the hole and serial number (as depicted in the image below)(if you are not familiar with how to safely complete this process, simply return the unit to Chart).

*Note:* If you have any question whether your unit is subject to this recall, please contact us immediately via email at MVEtankservice@chartindustries.com.



#### **INSTRUCTIONS FOR RETURN:**

Please identify the serial number for all units you possess that are subject to the recall, and identify the location of said products. Chart will arrange for return shipping. For processing the return, please send an email to: MVEtankservice@chartindustries.com.

If you elect to drill a hole in a retired unit, please email photograph(s) of the drilled unit(s) showing the hole and serial number(s).

If you no longer possess affected units, please immediately forward this communication to any entity in possession of affected units. Please also

**immediately** provide Chart with the location and contact information of the person/entity currently in possession of the unit(s). As soon as Chart has this information, Chart will coordinate replacement of the unit with you.

Chart/MVE is proud to be a trusted partner for our customers since 1963. We have worked hard throughout this time to set and uphold the quality standards in the industry. Please accept our sincerest apologies for any inconvenience that replacing these units may cause. We ask for your support as we work hand in hand with you to expeditiously remove the impacted units from service and restore your operations with quality units and trust in our products and processes.

Sincerely,

Earl Lawson

President, Chart Biomedical

Eare Causer

2200 Airport Industrial Drive, Suite 500

Ball Ground, GA 30107

From within the U.S. Dial 1-800-482-2473

From outside the U.S. dial +1-770-225-4428

## **ALUMINUM DEWARS RECALLED**

# **Serial Numbers:**

NPB201750xxxxx NPB201751xxxxx NPB201752xxxxx NPB201801xxxxx NPB201802xxxxx NPB201803xxxxx NPB201804xxxxx NPB201805xxxxx NPB201806xxxxx NPB201807xxxxx NPB201809xxxxx NPB201809xxxxx

## Part Number Description

10508967	FNL CRYSHPR W/PPSC
10647692	FNL CRYO-MOO-VER 3100
10650197	FNL CRYO-SYSTEM 4000
10650200	FNL CRYO-SYSTEM 2000
10718067	FNL CRYO-SYSTEM 6000
10719924	FNL XC47/11-6 W/11" CANISTERS
10725435	FNL XC47/11-10
10740281	FNL LAB 10
10742999	FNL BIO CANE 34
10743027	FNL XC34/18
10743035	FNL SC33/26
10743351	FNL BIO CANE 47
10770455	FNL LOCATOR 8 +
10770463	FNL LOCATOR 4 +
10770498	FNL LOCATOR JR +
10777411	FNL CRYSHPR XC IATA W/PPSC
10801514	FNL LOCATOR 6 +
10817330	FNL SC4/3V
10829621	FNL SC36/32
10979435	FNL CRYOSHIPPER MINI 4" OPEN
11005317	FNL SC20/12V
11015195	FNL CRYSHPR XC W/PPSC
11016649	FNL XC20/3V
11025650	FNL B20/20 BRYMILL
11037378	FNL CRYOSHIPPER MINI 2" OPEN
11037386	FNL CRYOSHIPPER MINI NO ABS
11064641	FNL CRYO-SYSTEM 2000 W/O C/C
11064659	FNL CRYO-SYSTEM 4000 W/O C/C
11064667	FNL CRYO-SYSTEM 6000 W/O C/C
11366673	FNL CRYOSHIPPER 2000

11378746 FNL XC47/11-6 W/O CANISTERS 11379221 FNL SC4/3V NO LABELS FNL DOBLE 20 SERIES 2001 11492151 11497948 FNL DOBLE 34 SERIES 2001 FNL DOBLE 47/6 SERIES 2001 11498684 11507886 FNL DOBLE 11 SERIES 2001 11527730 FNL DOBLE 28XC 11531544 FNL LAB 4 NO LABELS 11531552 FNL LAB 10 NO LABELS 11531579 FNL LAB 30 NO LABELS FNL CRYO-SYSTEM 2000 NO LABELS 11531587 FNL CRYO-SYSTEM 4000 NO LABELS 11531595 11531608 FNL CRYO-SYSTEM 6000 NO LABELS 11536492 FNL XC34/18 ABS 11541339 FNL LAB 50 NO LABELS 11549357 FNL MIL XC20 NO LABELS 11554308 FNL DOBLE 47-10 11647424 FNL LAB 5 NO LABELS 11736179 FNL SPECTRUM 8 FNL SPECTRUM 12 11736187 FNL SPECTRUM 20 W/11" CNSTR 11736208 FNL SPECTRUM 35 W/11" CNSTR 11736216 11778168 FNL ARCTIC EXPRESS DUAL 10 FNL ARCTIC EXPRESS DUAL 19 11778176 11778184 FNL ARCTIC EXPRESS DUAL 28 11778230 FNL ARCTIC EXPRESS IATA 11838853 FNL SPECTRUM 3 FNL CRYO-SYSTEM 750 11886450 11901162 FNL MILLENNIUM II W/11" CNSTR 11901584 FNL SC11/7 W/135MM CANS 11913497 FNL LAB 47 FNL CRYSHPR CRYOLIFE 11929481 13101839 FNL XC34/18 PLUS W/WELD CNSTRS FNL SPECTRUM 8 NO CANISTERS 13112765 FNL XC47/11-6 W/XL SQ CNSTRS 13346449 13490281 CFG SPECTRUM 35 LITER EUROPE 13492631 FNL LAB 20 13512890 FNL LAB 20 NO LABELS 13571523 FNL ET40-10 **CFG 47 LITER EUROPE** 13635355 13798111 FNL THERMO 20 14656887 FNL LOCATOR 6 FNL BIO CANE 73 14670312 FNL ARCTIC EXPRESS 20 14670321 14796775 FNL CRYOSYSTEMS 6000 FULL AUTO 14893049 FNL XC-20 SELECT SIRES 20638227 FNL SEMEX-20 FNL CRI-20 20639136 20640380 FNL ALTA-20 FNL MVE-20 20640771 20668449 FNL SPECTRUM 20 NO CANISTERS 20677294 FNL CRYSHPR XC IATA W/O SC

FNL DOBLE 47/6 CEVA

20726502

20726517 FNL DOBLE 11 CEVA 20789363 FNL CRYOSYSTEMS 6000 AF-GB-MDD 20861774 FNL SC20/20 SIG 20875477 FNL ARCTIC EXPRESS 20 DL FNL CRYSHPR W/PPSC W/ SQ RACK 20925284 20937586 FNL DOBLE 22 20947690 FNL XC20 SIG ABS CFG XC20 SIG 20997044 20997986 FNL XC20 SIG 20999720 FNL XC20 SIG SELECT SIRES FNL SPECTRUM 35 W/O CNSTR 21007470 FNL MVE WS-20 SIG 21010336 21036452 FNL XC20 SIG ALTA 21036453 FNL XC20 SIG GENEX 21036454 FNL XC20 SIG SEMEX 21099759 FNL SEXING TECH- 20 21115938 FNL BIO CANE 20 SIG 21134753 FNL SC4/3V W/PPSC 21134754 FNL SC4/2V W/PPSC 21178824 FNL SC20/20 SIG PREDICTIVE BIO 9917989 FNL SC8/5 9918009 FNL SC3/3 9918039 FNL XC43/28 FNL XC33/22 9918069 9918079 FNL LAB 5 9918099 FNL LAB 30 9918109 FNL LAB 50 9918449 FNL XC32/8 9918479 FNL SC4/2V 9918499 FNL SC11/7 9918539 FNL XC22/5 9918999 FNL XC47/11 MICRO-MERETICS FNL XC21/6E 9920579 9920759 FNL THERMO 5 9920769 FNL THERMO 10 9920789 FNL THERMO 30 9920879 **FNL LOCATOR 8** 9920899 **FNL LOCATOR 4** 9920919 FNL LOCATOR JR 9921849 FNL ARCTIC EXPRESS 10 9921979 FNL ARCTIC EXPRESS 5 9922219 FNL LAB 4 9922339 FNL SC2/1V