



## URGENT MEDICAL DEVICE FIELD SAFETY CORRECTIVE ACTION

### Hookup Accessories used with MEDIVATORS ADVANTAGE PLUS, ADVANTAGE PLUS Pass Thru, and ISA Automated Endoscope Reprocessors

November 20, 2023

**ATTN: GI DIRECTOR OR STERILE PROCESSING DEPARTMENT**

Dear Valued STERIS Customer:

STERIS is voluntarily conducting a recall for hookup accessories used with the MEDIVATORS™ ADVANTAGE PLUS™, ADVANTAGE PLUS Pass Thru, and ISA Automated Endoscope Reprocessors **distributed between July 21 – November 13, 2023**. Our records indicate that your facility purchased one or more of the hookups impacted by this recall. Please reference Attachment A to this letter for the full list of affected models and UDIs.

**Description of the product** – Hookups are intended to provide the interface between an Automated Endoscope Reprocessor (AER) and endoscopes, which allow perfusion of liquids dispensed by the AER to contact and high-level disinfect the channels of flexible endoscopes.

**Description of the problem** – Through receipt of Customer complaints, it has been determined that the identification label on affected hookups could become illegible over time due to the label's ink coming off. If the information on the hookup label were to become illegible, this may result in an inconvenience to users, or in some cases where a user facility does not have a backup hookup, potentially prevent the user from being able to effectively identify the hookup (i.e., delineate one affected hookup from another). This may result in the user electing to not reprocess certain endoscopes, which could result in postponement of patient procedures.

**STERIS Action** – STERIS will provide either a replacement hookup or replacement label for all affected hookups.

**User Action** – Please ensure the following steps are completed:

1. Please inspect your on-hand inventory for product affected by this recall. You should not attempt to use a hookup that cannot be effectively identified.
2. Please complete the Medical Device Recall Response Form included with this letter, and return the completed Response Form via email to: [Regulatory\\_Compliance@STERIS.com](mailto:Regulatory_Compliance@STERIS.com) or via fax to 440-392-8963. **IMPORTANT NOTE:** No replacement hookup(s) or label(s) will be provided until a completed Response Form is received.

We apologize for any inconvenience this matter may cause, and as always, STERIS is dedicated to supporting our products and valued Customers. If you have questions regarding this matter, please contact [Regulatory\\_Compliance@STERIS.com](mailto:Regulatory_Compliance@STERIS.com).

Sincerely,

A handwritten signature in blue ink that reads "Michelle LaVan".

Michelle LaVan  
Lead, Quality & Regulatory Compliance Specialist  
STERIS

## MEDICAL DEVICE RECALL RESPONSE FORM

### **RESPONSE IS REQUIRED**

Facility Name: \_\_\_\_\_

Street Address: \_\_\_\_\_

City, State, Zip/Postal Code: \_\_\_\_\_

### **Hookup Accessories used with ADVANTAGE PLUS, ADVANTAGE PLUS Pass Thru, and ISA Automated Endoscope Reprocessors**

IMPORTANT: Replacement hookup(s) or replacement label(s) will be provided based on affected product shipped to your facility. Please fill out this form in its entirety and return to STERIS to initiate the replacement process.

\_\_\_\_\_  
Printed Name and Title

\_\_\_\_\_  
\*Signature and Date

\*By signing this Response Form, you are confirming that your facility will not attempt to use a hookup that cannot be effectively identified.

**Please fill out in entirety, scan, and return via email to  
Regulatory\_Compliance@STERIS.com or via fax to 440-392-8963.**

## Attachment A – List of affected Hookup models and UDIs

| Model # | Model Description         | UDI            |
|---------|---------------------------|----------------|
| 2-8-001 | HU ZERO CHANNEL NO LT     | 00677964068515 |
| 2-8-010 | ASSY, HOOK-UP, OLYMPUS    | 00677964068478 |
| 2-8-014 | HOOKUP, OPTIM LT          | 00677964095177 |
| 2-8-030 | ASSY, HOOKUP, PENTAX-0K   | 00677964071348 |
| 2-8-040 | ASSY, HKUP STORZ 0K       | 00677964072284 |
| 2-8-070 | ASSY, HU VISION SCIENCE   | 00677964094637 |
| 2-8-112 | ASSY, HOOKUP OLMPUS 1K    | 00677964071300 |
| 2-8-120 | ASSY, HOOKUP FUJI 1K      | 00677964102080 |
| 2-8-210 | ASSY, HOOK-UP OLYMPUS 2K  | 00677964068485 |
| 2-8-213 | ASSY, HOOKUP, OLYMPUS-2K  | 00677964070648 |
| 2-8-218 | ASSY HOOKUP OLY 2K        | 00677964072369 |
| 2-8-236 | ASSY, HOOKUP FUJINON 2K   | 00677964072376 |
| 2-8-244 | ASSY, HOOKUP, STORZ 2K    | 00677964071836 |
| 2-8-245 | ASSY, HOOK-UP, STORZ 2K   | 00677964071195 |
| 2-8-246 | ASSY, HKUP, PENTAX-2K     | 00677964084058 |
| 2-8-260 | ASSY, HOOK-UP WOLF 2-K    | 00677964105951 |
| 2-8-280 | ASSY HOOKUP               | 00677964073069 |
| 2-8-312 | ASSY, HOOK-UP OLYMPUS-3K  | 00677964074400 |
| 2-8-330 | ASSY, HOOK-UP, PENTAX-3K  | 00677964073816 |
| 2-8-332 | ASSEMBLY, HOOK-UP, PENTAX | 00677964091483 |
| 2-8-334 | ASSEMBLY, HOOK-UP, PENTAX | 00677964073830 |
| 2-8-335 | ASSY, HOOKUP, OLYMPUS, 3K | 00677964071263 |
| 2-8-336 | ASSY, HKUP FUJI 70 BRONC  | 00677964091506 |
| 2-8-338 | ASSY, HU, FUJI BRONCH U/S | 00677964074141 |
| 2-8-421 | ASSY, HOOKUP, FUJIFILM    | 00677964071270 |
| 2-8-428 | ASSY, ED-580XX            | 00677964104923 |
| 2-8-431 | ASSY, HOOKUP, PENTAX      | 00677964071706 |
| 2-8-432 | ASSY, HOOK-UP, PENTAX     | 00677964072239 |
| 2-8-460 | HOOKUP, ASSY WOLF 4K      | 00677964092756 |
| 2-8-511 | ASSY, HOOK-UP, OLY        | 00677964066726 |
| 2-8-513 | ASSY HOOKUP OLY 4K ESOPHA | 00677964102394 |
| 2-8-522 | ASSY, FUJINON 4K-JET      | 40677964025387 |
| 2-8-523 | ASSY, HKUP, FUJI 4K SER 5 | 00677964074158 |
| 2-8-526 | ASSY, HOOKUP FUJI 4K+JET  | 00677964071904 |
| 2-8-527 | ASSY, FUJI 700 NO FWJ     | 00677964111648 |
| 2-8-528 | HOOK, EI-580BT            | 40677964025394 |
| 2-8-534 | ASSY, HOOK-UP, PEN-4K+JET | 00677964073632 |
| 2-8-537 | ASSY, HKUP, PTX90 4K+JET  | 00677964068492 |
| 2-8-539 | HOOK, PTX I20 FWJ         | 40677964039827 |
| 2-8-540 | ASSY-HOOK-UP, STORZ, 5K   | 00677964091223 |
| 2-8-611 | ASSY, HOOK-UP, OLY        | 00677964066719 |
| 2-8-612 | ASSY, HOOKUP, PENTAX      | 00677964073571 |
| 2-8-614 | ASM HU OLY TGF-UC180J     | 00677964101984 |
| 2-8-617 | ASSY, HKUP, OLY 2 BIOPSY  | 00677964073427 |
| 2-8-623 | ASSY, HKUP, FUJI-5 4K+JET | 00677964091476 |
| 2-8-625 | ASSY, HOOKUP, EG-580UT    | 00677964077678 |
| 2-8-626 | ASSY, HOOKUP, EG-580UR    | 00677964077685 |
| 2-8-627 | ASSY, HKUP, FUJI 700 SER  | 00677964083006 |
| 2-8-634 | ASSY, HOOKUP PENTAX DUAL  | 00677964074189 |
| 2-8-639 | HOOK, EG3X-J10            | 40677964006560 |
| 2-8-640 | ASSY, HOOK-UP STORZ 6K    | 00677964074417 |

| <b>Model #</b> | <b>Model Description</b>  | <b>UDI</b>     |
|----------------|---------------------------|----------------|
| 2-8-711        | ASSY, HOOKUP, OLYMPUS 5K  | 00677964070662 |
| 2-8-718        | ASSY, HOOKUP, OLY 160/180 | 00677964070655 |
| 2-8-740        | ASSY, HOOK-UP, STORZ      | 00677964102479 |
| 2-8-811        | ASSY HKUP OLY-ECHO+LIFT   | 00677964070624 |
| 2-8-901        | KIT TEE/TOE/RUS HU ADV+   | 00677964046926 |
| 5-8-010        | ASSY, HOOKUP DBL LT       | 00677964071874 |
| 5-8-030        | ASSY, HOOKUP              | 00677964094644 |
| 5-8-040        | ASSY HOOKUP STORZ 0K DBL  | 00677964072772 |
| 5-8-050        | ASSY, HOOKUP, DILATOR     | 00677964072963 |
| 5-8-112        | HOOKUP 2-SCOPE 1-K        | 00677964072307 |
| 5-8-190        | HOOK DBL, IS, ION         | 40677964025493 |
| 5-8-210        | ASSY, HOOKUP              | 00677964080333 |
| 5-8-245        | HOOKUP KARLSTORZ-2K       | 00677964102493 |
| 2-8-001HAN     | HU ZERO CHANNEL NO LT     | 00677964064739 |
| 2-8-002HAN     | ADVANTAGE TEST BLOCK      | 00677964046568 |
| 2-8-010HAN     | ASSEMBLY, HOOK-UP,OLYMPUS | 00677964001802 |
| 2-8-030HAN     | ASSY,HOOKUP,PENTAX-OK     | 00677964037399 |
| 2-8-112HAN     | ASSY, HOOK-UP OLYMPUS 1K  | 00677964001819 |
| 2-8-160HAN     | ASSY, ADV HOOKUP, WOLF    | 00677964037443 |
| 2-8-210HAN     | ASSY, HOOK-UP OLYMPUS 2K  | 00677964001826 |
| 2-8-213HAN     | ASSY, HOOKUP, OLYMPUS-2K  | 00677964001871 |
| 2-8-218HAN     | ASSY HOOKUP OLY 2K        | 00677964061707 |
| 2-8-245HAN     | ASSY, HOOKUP, STORZ-2K    | 00677964037511 |
| 2-8-330HAN     | ASSY, HOOK-UP, PENTAX-3K  | 00677964001970 |
| 2-8-334HAN     | ASSEMBLY, HOOK-UP, PENTAX | 00677964001994 |
| 2-8-335HAN     | ASSY,HOOKUP,OLYMPUS       | 00677964037634 |
| 2-8-336HAN     | ASSY,HKUP FUJI 70 BRONC   | 00677964045387 |
| 2-8-338HAN     | ASSY, HU, FUJI BRONCH U/S | 00677964053443 |
| 2-8-421HAN     | ASSY, HOOKUP, FUJINON     | 00677964037658 |
| 2-8-428CAS     | HOOK, CAS, ED-580XX       | 40677964024663 |
| 2-8-428HAN     | HOOK, HAN, ED-580XX       | 40677964024670 |
| 2-8-431CAS     | ASSY, CASSETTE, PENTAX    | 00677964002045 |
| 2-8-431HAN     | ASSY, HOOK-UP, PENTAX     | 00677964002052 |
| 2-8-511HAN     | ASSY, HOOK-UP, OLY        | 00677964002144 |
| 2-8-526HAN     | ASSY, HOOKUP FUJI 4K+JET  | 00677964053764 |
| 2-8-536HAN     | HOOKUP,PENTAX,70ULT, 5K   | 00677964044168 |
| 2-8-537CAS     | ASSY,CASS,PTX 90 4K+JET   | 00677964037801 |
| 2-8-537HAN     | ASSY HKUP,PTX90 4K+JET    | 00677964037818 |
| 2-8-611HAN     | ASSY, HOOK-UP, OLY        | 00677964002236 |
| 2-8-626CAS     | ASSY, CASS, FUJINON 580UR | 00677964077661 |
| 2-8-627HAN     | ASSY, HAN, FUJI 700 SER   | 00677964089701 |
| 2-8-634CAS     | ASSY,CASS PENTAX DUAL CH  | 00677964054815 |
| 2-8-634HAN     | ASSY,HOOKUP PENTAX DUAL   | 00677964054808 |
| 2-8-640HAN     | ASSY, HOOK-UP STORZ 6K    | 00677964087097 |
| 2-8-711HAN     | ASSY, HOOKUP, OLYMPUS-5K  | 00677964037924 |
| 2-8-718HAN     | ASSY,HOOK-UP, OLY 160/180 | 00677964064777 |
| 2-8-811HAN     | ASSY HKUP OLY-ECHO W/LIFT | 00677964053542 |
| 2-8-902H       | HU BLOCK ULTRASOUND PROBE | 00677964061509 |
| 5-8-030HAN     | ASSY, HOOKUP              | 00677964046681 |
| 5-8-050HAN     | ASSY,HOOKUP, DILATOR      | 00677964040306 |
| 5-8-112HAN     | ASSY, HOOKUP 2-SCOPE 1-K  | 00677964051548 |
| 5-8-244HAN     | ASSY, HKUP 2K DBL         | 00677964045165 |



## URGENT MEDICAL DEVICE FIELD SAFETY CORRECTIVE ACTION

### Hookup Accessories used with MEDIVATORS ADVANTAGE PLUS, ADVANTAGE PLUS Pass Thru, and ISA Automated Endoscope Reprocessors

November 20, 2023

**ATTN: PURCHASING DEPARTMENT**

Dear Valued STERIS Distributor:

STERIS is voluntarily conducting a recall for hookup accessories used with the MEDIVATORS™ ADVANTAGE PLUS™, ADVANTAGE PLUS Pass Thru, and ISA Automated Endoscope Reprocessors **distributed between July 21 – November 13, 2023**. Our records indicate that your facility purchased one or more of the hookups impacted by this recall. Please reference Attachment A to this letter for the full list of affected models and UDIs.

**Description of the product** – Hookups are intended to provide the interface between an Automated Endoscope Reprocessor (AER) and endoscopes, which allow perfusion of liquids dispensed by the AER to contact and high-level disinfect the channels of flexible endoscopes.

**Description of the problem** – Through receipt of Customer complaints, it has been determined that the identification label on affected hookups could become illegible over time due to the label's ink coming off. If the information on the hookup label were to become illegible, this may result in an inconvenience to users, or in some cases where a user facility does not have a backup hookup, potentially prevent the user from being able to effectively identify the hookup (i.e., delineate one affected hookup from another). This may result in the user electing to not reprocess certain endoscopes, which could result in postponement of patient procedures.

**STERIS Action** – STERIS will provide either a replacement hookup or replacement label for all affected hookups.

**User Action** – Please ensure the following steps are completed:

1. Please inspect your on-hand inventory for product affected by this recall. You should not attempt to use a hookup that cannot be effectively identified.
2. Please complete the Medical Device Recall Response Form included with this letter, and return the completed Response Form via email to: [Regulatory\\_Compliance@STERIS.com](mailto:Regulatory_Compliance@STERIS.com) or via fax to 440-392-8963. **IMPORTANT NOTE:** No replacement hookup(s) or label(s) will be provided until a completed Response Form is received.
3. If you have further distributed this product, please identify your Customer(s), and notify them of this field recall.

We apologize for any inconvenience this matter may cause, and as always, STERIS is dedicated to supporting our products and valued Customers. If you have questions regarding this matter, please contact [Regulatory\\_Compliance@STERIS.com](mailto:Regulatory_Compliance@STERIS.com).

Sincerely,

A handwritten signature in blue ink that reads "Michelle LaVan".

Michelle LaVan  
Lead, Quality & Regulatory Compliance Specialist  
STERIS

## MEDICAL DEVICE RECALL RESPONSE FORM

### **RESPONSE IS REQUIRED**

Facility Name: \_\_\_\_\_

Street Address: \_\_\_\_\_

City, State, Zip/Postal Code: \_\_\_\_\_

### **Hookup Accessories used with ADVANTAGE PLUS, ADVANTAGE PLUS Pass Thru, and ISA Automated Endoscope Reprocessors**

IMPORTANT: Replacement hookup(s) or replacement label(s) will be provided based on affected product shipped to your facility. Please fill out this form in its entirety and return to STERIS to initiate the replacement process.

\_\_\_\_\_  
Printed Name and Title

\_\_\_\_\_  
\*Signature and Date

\*By signing this Response Form, you are confirming that your facility will not attempt to use a hookup that cannot be effectively identified.

**Please fill out in entirety, scan, and return via email to  
Regulatory\_Compliance@STERIS.com or via fax to 440-392-8963.**

## Attachment A – List of affected Hookup models and UDIs

| Model # | Model Description         | UDI            |
|---------|---------------------------|----------------|
| 2-8-001 | HU ZERO CHANNEL NO LT     | 00677964068515 |
| 2-8-010 | ASSY, HOOK-UP, OLYMPUS    | 00677964068478 |
| 2-8-014 | HOOKUP, OPTIM LT          | 00677964095177 |
| 2-8-030 | ASSY, HOOKUP, PENTAX-0K   | 00677964071348 |
| 2-8-040 | ASSY, HKUP STORZ 0K       | 00677964072284 |
| 2-8-070 | ASSY, HU VISION SCIENCE   | 00677964094637 |
| 2-8-112 | ASSY, HOOKUP OLMPUS 1K    | 00677964071300 |
| 2-8-120 | ASSY, HOOKUP FUJI 1K      | 00677964102080 |
| 2-8-210 | ASSY, HOOK-UP OLYMPUS 2K  | 00677964068485 |
| 2-8-213 | ASSY, HOOKUP, OLYMPUS-2K  | 00677964070648 |
| 2-8-218 | ASSY HOOKUP OLY 2K        | 00677964072369 |
| 2-8-236 | ASSY, HOOKUP FUJINON 2K   | 00677964072376 |
| 2-8-244 | ASSY, HOOKUP, STORZ 2K    | 00677964071836 |
| 2-8-245 | ASSY, HOOK-UP, STORZ 2K   | 00677964071195 |
| 2-8-246 | ASSY, HKUP, PENTAX-2K     | 00677964084058 |
| 2-8-260 | ASSY, HOOK-UP WOLF 2-K    | 00677964105951 |
| 2-8-280 | ASSY HOOKUP               | 00677964073069 |
| 2-8-312 | ASSY, HOOK-UP OLYMPUS-3K  | 00677964074400 |
| 2-8-330 | ASSY, HOOK-UP, PENTAX-3K  | 00677964073816 |
| 2-8-332 | ASSEMBLY, HOOK-UP, PENTAX | 00677964091483 |
| 2-8-334 | ASSEMBLY, HOOK-UP, PENTAX | 00677964073830 |
| 2-8-335 | ASSY, HOOKUP, OLYMPUS, 3K | 00677964071263 |
| 2-8-336 | ASSY, HKUP FUJI 70 BRONC  | 00677964091506 |
| 2-8-338 | ASSY, HU, FUJI BRONCH U/S | 00677964074141 |
| 2-8-421 | ASSY, HOOKUP, FUJIFILM    | 00677964071270 |
| 2-8-428 | ASSY, ED-580XX            | 00677964104923 |
| 2-8-431 | ASSY, HOOKUP, PENTAX      | 00677964071706 |
| 2-8-432 | ASSY, HOOK-UP, PENTAX     | 00677964072239 |
| 2-8-460 | HOOKUP, ASSY WOLF 4K      | 00677964092756 |
| 2-8-511 | ASSY, HOOK-UP, OLY        | 00677964066726 |
| 2-8-513 | ASSY HOOKUP OLY 4K ESOPHA | 00677964102394 |
| 2-8-522 | ASSY, FUJINON 4K-JET      | 40677964025387 |
| 2-8-523 | ASSY, HKUP, FUJI 4K SER 5 | 00677964074158 |
| 2-8-526 | ASSY, HOOKUP FUJI 4K+JET  | 00677964071904 |
| 2-8-527 | ASSY, FUJI 700 NO FWJ     | 00677964111648 |
| 2-8-528 | HOOK, EI-580BT            | 40677964025394 |
| 2-8-534 | ASSY, HOOK-UP, PEN-4K+JET | 00677964073632 |
| 2-8-537 | ASSY, HKUP, PTX90 4K+JET  | 00677964068492 |
| 2-8-539 | HOOK, PTX I20 FWJ         | 40677964039827 |
| 2-8-540 | ASSY-HOOK-UP,STORZ,5K     | 00677964091223 |
| 2-8-611 | ASSY, HOOK-UP, OLY        | 00677964066719 |
| 2-8-612 | ASSY, HOOKUP, PENTAX      | 00677964073571 |
| 2-8-614 | ASM HU OLY TGF-UC180J     | 00677964101984 |
| 2-8-617 | ASSY, HKUP, OLY 2 BIOPSY  | 00677964073427 |
| 2-8-623 | ASSY, HKUP, FUJI-5 4K+JET | 00677964091476 |
| 2-8-625 | ASSY, HOOKUP, EG-580UT    | 00677964077678 |
| 2-8-626 | ASSY, HOOKUP, EG-580UR    | 00677964077685 |
| 2-8-627 | ASSY, HKUP, FUJI 700 SER  | 00677964083006 |
| 2-8-634 | ASSY, HOOKUP PENTAX DUAL  | 00677964074189 |
| 2-8-639 | HOOK, EG3X-J10            | 40677964006560 |
| 2-8-640 | ASSY, HOOK-UP STORZ 6K    | 00677964074417 |

| <b>Model #</b> | <b>Model Description</b>  | <b>UDI</b>     |
|----------------|---------------------------|----------------|
| 2-8-711        | ASSY, HOOKUP, OLYMPUS 5K  | 00677964070662 |
| 2-8-718        | ASSY, HOOKUP, OLY 160/180 | 00677964070655 |
| 2-8-740        | ASSY, HOOK-UP, STORZ      | 00677964102479 |
| 2-8-811        | ASSY HKUP OLY-ECHO+LIFT   | 00677964070624 |
| 2-8-901        | KIT TEE/TOE/RUS HU ADV+   | 00677964046926 |
| 5-8-010        | ASSY, HOOKUP DBL LT       | 00677964071874 |
| 5-8-030        | ASSY, HOOKUP              | 00677964094644 |
| 5-8-040        | ASSY HOOKUP STORZ 0K DBL  | 00677964072772 |
| 5-8-050        | ASSY, HOOKUP, DILATOR     | 00677964072963 |
| 5-8-112        | HOOKUP 2-SCOPE 1-K        | 00677964072307 |
| 5-8-190        | HOOK DBL, IS, ION         | 40677964025493 |
| 5-8-210        | ASSY, HOOKUP              | 00677964080333 |
| 5-8-245        | HOOKUP KARLSTORZ-2K       | 00677964102493 |
| 2-8-001HAN     | HU ZERO CHANNEL NO LT     | 00677964064739 |
| 2-8-002HAN     | ADVANTAGE TEST BLOCK      | 00677964046568 |
| 2-8-010HAN     | ASSEMBLY, HOOK-UP,OLYMPUS | 00677964001802 |
| 2-8-030HAN     | ASSY,HOOKUP,PENTAX-OK     | 00677964037399 |
| 2-8-112HAN     | ASSY, HOOK-UP OLYMPUS 1K  | 00677964001819 |
| 2-8-160HAN     | ASSY, ADV HOOKUP, WOLF    | 00677964037443 |
| 2-8-210HAN     | ASSY, HOOK-UP OLYMPUS 2K  | 00677964001826 |
| 2-8-213HAN     | ASSY, HOOKUP, OLYMPUS-2K  | 00677964001871 |
| 2-8-218HAN     | ASSY HOOKUP OLY 2K        | 00677964061707 |
| 2-8-245HAN     | ASSY, HOOKUP, STORZ-2K    | 00677964037511 |
| 2-8-330HAN     | ASSY, HOOK-UP, PENTAX-3K  | 00677964001970 |
| 2-8-334HAN     | ASSEMBLY, HOOK-UP, PENTAX | 00677964001994 |
| 2-8-335HAN     | ASSY,HOOKUP,OLYMPUS       | 00677964037634 |
| 2-8-336HAN     | ASSY,HKUP FUJI 70 BRONC   | 00677964045387 |
| 2-8-338HAN     | ASSY, HU, FUJI BRONCH U/S | 00677964053443 |
| 2-8-421HAN     | ASSY, HOOKUP, FUJINON     | 00677964037658 |
| 2-8-428CAS     | HOOK, CAS, ED-580XX       | 40677964024663 |
| 2-8-428HAN     | HOOK, HAN, ED-580XX       | 40677964024670 |
| 2-8-431CAS     | ASSY, CASSETTE, PENTAX    | 00677964002045 |
| 2-8-431HAN     | ASSY, HOOK-UP, PENTAX     | 00677964002052 |
| 2-8-511HAN     | ASSY, HOOK-UP, OLY        | 00677964002144 |
| 2-8-526HAN     | ASSY, HOOKUP FUJI 4K+JET  | 00677964053764 |
| 2-8-536HAN     | HOOKUP,PENTAX,70ULT, 5K   | 00677964044168 |
| 2-8-537CAS     | ASSY,CASS,PTX 90 4K+JET   | 00677964037801 |
| 2-8-537HAN     | ASSY HKUP,PTX90 4K+JET    | 00677964037818 |
| 2-8-611HAN     | ASSY, HOOK-UP, OLY        | 00677964002236 |
| 2-8-626CAS     | ASSY, CASS, FUJINON 580UR | 00677964077661 |
| 2-8-627HAN     | ASSY, HAN, FUJI 700 SER   | 00677964089701 |
| 2-8-634CAS     | ASSY,CASS PENTAX DUAL CH  | 00677964054815 |
| 2-8-634HAN     | ASSY,HOOKUP PENTAX DUAL   | 00677964054808 |
| 2-8-640HAN     | ASSY, HOOK-UP STORZ 6K    | 00677964087097 |
| 2-8-711HAN     | ASSY, HOOKUP, OLYMPUS-5K  | 00677964037924 |
| 2-8-718HAN     | ASSY,HOOK-UP, OLY 160/180 | 00677964064777 |
| 2-8-811HAN     | ASSY HKUP OLY-ECHO W/LIFT | 00677964053542 |
| 2-8-902H       | HU BLOCK ULTRASOUND PROBE | 00677964061509 |
| 5-8-030HAN     | ASSY, HOOKUP              | 00677964046681 |
| 5-8-050HAN     | ASSY,HOOKUP, DILATOR      | 00677964040306 |
| 5-8-112HAN     | ASSY, HOOKUP 2-SCOPE 1-K  | 00677964051548 |
| 5-8-244HAN     | ASSY, HKUP 2K DBL         | 00677964045165 |