

## **URGENT FIELD SAFETY NOTICE**

## **RECALL of CRP FS**

## Deviation in recovery in lower concentration range

Date 2023-06-26

Product / Lot CRP FS Cat. No.: 1 7002 ...

Lot: 33801/50293937, Kit Lot: 60158540, 60158618,

60158731, 60158786, 60158922, 60158937, 60158981, 60159004, 60159088, 60159295

Lot: 34108, Kit Lot: 60158891, 60159263, 60159308,

60159309, 60159320, 60159324, 60159498, 60159595, 60159619, 60159634, 60159840,

60159841, 60159892, 60160225

Lot: 34574, Kit Lot: 60159445

Explanation Based on the analysis of a customer complaint the a.m. lots of CRP FS

> have been reviewed. The internal inspection confirms incorrect recovery of low concentrated samples (<10 mg/L), which leads to false high results. The cause has been traced back to a specific raw material.

Impact on patient

results

CRP FS measurements of patient samples lead to deviating results.

The a.m. reagent lots are withdrawn from the market. Measures

> Please immediately inform all involved users about this recall. Immediately stop further measurements and replace the reagent lots at the customer's site. Existing stock must be disposed of.

DiaSys offers replacement for the affected kits.

Determinations performed using a.m. lots of CRP FS shall be assessed

with care and repeated if implausible results occur.

DiaSys has announced this field safety notice to the relevant authorities of the European Union. Customers outside the EU are asked to handle necessary announcements to authorities in their countries.

Under current regulations, we are obliged to provide a complete chain of evidence of all corrective measures for our products. For this reason, we ask you to fill in and sign the attached confirmation of receipt as proof that you have received and communicated this information to all concerned customers. Please return it by fax or as scan until 2023-07-31.





## **URGENT FIELD SAFETY NOTICE**

Please accept our sincere apologies for the inconvenience caused. In case you have any questions, please do not hesitate to contact us.

Best regards,

Tim Fleissner

Strategic Product Management Reagents

